

# Supercritical Fluids in Reaction Engineering

## Syringe Pump Application Note AN5

### Introduction

The benefits of Supercritical Fluids in Reaction Engineering applications can greatly exceed the initial equipment outlay. Research in reactions involving supercritical fluids, e.g., SFCO<sub>2</sub>, has shown it possible to obtain the following advantages:


- better product uniformity
- faster reaction rates
- improved selectivity
- greater energy savings compared to evaporation of traditional bulk solvents
- more environmentally benign
- non-flammable
- non-toxic

Teledyne Isco syringe pumps provide accurate pulse-less flow and are excellent CO<sub>2</sub> pumps or Supercritical Fluid pumps. They are commonly used in SF reaction engineering research at lab and pilot scales.

### Theory

Supercritical fluids are very dense gases with many properties superior to liquids or solvents. While there are many fluids that can be used in their supercritical state, CO<sub>2</sub> is the one most often used because it is considered environmentally friendly, and its critical temperature and operating pressures are relatively easy to work with. In the supercritical state, molecular forces that give liquids their particular properties of surface tension, viscosity or slower diffusion are altered. Molecules do not “stick together” as well, so viscosities are lower with higher diffusion rates. Such enhanced properties can be beneficial to the reaction process since mixing is improved thereby improving distribution and enhancing product quality.

A phase diagram for CO<sub>2</sub>, shown in Figure 1, displays the relationship between pressure and temperature. When the conditions of pressure and temperature are altered, the phases of CO<sub>2</sub> can be changed to a solid, liquid, or gas. However, when above the critical temperature, TC, CO<sub>2</sub> becomes supercritical, and can no longer be changed back into liquid by increasing the pressure. In this state, CO<sub>2</sub> will remain a gas-like fluid even though it may be approaching the density of a liquid at very high pressures. Its supercritical properties include solvating power similar to liquids, the penetrating or diffusion properties of a gas, and a “zero” surface tension.



**Figure 1: Phase Diagram of Carbon Dioxide**

CO<sub>2</sub> is a naturally occurring component of our atmosphere. CO<sub>2</sub> is non-toxic in small amounts (consider your own breath), and not a volatile organic compound (VOC), hence not contributing to smog formation. CO<sub>2</sub> is also non-flammable—a great advantage over many conventional liquid solvents.

Supercritical CO<sub>2</sub> has a low viscosity and high diffusivity compared to the usual liquid solvents used as a medium for reactions. Low viscosity and high diffusivity cause the reagents (or soluble catalysts) to rapidly travel to all locations inside a reactor. With uniform conditions throughout the reactor, the reaction process can obtain desirable results: better product uniformity, faster reaction rates, and improved selectivity. It should be noted that, despite all the advantages, SFCO<sub>2</sub> is best at dissolving small, non-polar organic compounds, but has difficulty dissolving many polar or ionic compounds, or most large polymers (fluorinated oligomers are an exception). Solvating properties can be improved with the addition of small amounts of other fluids or modifiers. These can include additives such as surfactants, or modifiers such as ethanol or methanol.

After completing supercritical fluid reaction, SFCO<sub>2</sub> is depressurized to a gaseous state, and solutes generally precipitate and fall out of solution. That is, the reaction product(s) is no longer soluble in gaseous CO<sub>2</sub>. In the manufacture of chemical solids, often the product must be dried after the reaction process to remove solvents. A large amount of energy is used for drying in these conventional processes. With a reaction process involving SFCO<sub>2</sub>, the product is left dry after the reaction, and following depressurization; hence no further drying is required. This has also been shown in polymer manufacturing, where by using SFCO<sub>2</sub>, dry product can be obtained with a resultant significant energy savings

compared to the traditional drying procedures required to evaporate liquid solvents or water. Therefore, there can be very large savings in energy costs for a process utilizing CO<sub>2</sub>.

The CO<sub>2</sub> gas can then be recycled, another major cost reduction. When using CO<sub>2</sub> in a continuous manufacturing process, large amounts of CO<sub>2</sub>, which could function as greenhouse gas, are recycled instead of being released into the environment. Small CO<sub>2</sub> leaks do not harm the environment.

## Applications Using Supercritical CO<sub>2</sub>

**Catalysis** – Homogeneous catalysts are highly active and selective, while heterogeneous catalysts are less active but easier to separate and re-use. In 2006, Liotta, Eckert, Hallett, and Pollet reported techniques for recycling homogeneous catalysts using SFCO<sub>2</sub> by changing pressure of the CO<sub>2</sub> in the reaction. This allowed the homogeneity to be turned on and off. Another example comes from Toghiani et al, where SFCO<sub>2</sub> was used to oxidize unsaturated fatty acids to make diacids and epoxides. The reaction medium was SFCO<sub>2</sub>, which was completely oxidized.

**Nanotechnology** – SFCO<sub>2</sub> is becoming an enabling solvent for producing nanomaterials such as aerogels of Al<sub>2</sub>O<sub>3</sub>, SiO<sub>2</sub>, TiO<sub>2</sub>, and ZrO<sub>2</sub>. These nanomaterials, which often have new and exciting properties such as tunable pore sizes and high surface areas, have applications as biomaterials, and catalysts for fuel cells and solar cells. Also, SFCO<sub>2</sub> drying has been widely used to produce aerogels that exhibit a very high specific surface area and maintain the nanoarchitecture due to the zero surface tension.

**Pharmaceuticals and Biomedical Devices** – Supercritical fluids have emerged as the green solvents in the pharmaceutical industry for micronization of inhalable medicines, and separation of chiral enantiomers. SFCO<sub>2</sub> is also playing a role as a valuable tool in tissue engineering and preparing biomedical devices, as CO<sub>2</sub> is largely anti-bacterial and can be used in producing tissue scaffolds.

**Polymerization** – Previously, DeSimone and coworkers have shown that SFCO<sub>2</sub> is a promising alternative medium for free-radical, cationic, and step-growth polymerizations, and continuous processes. The free radical initiator AIBN was shown to have a higher efficiency in SFCO<sub>2</sub> than in benzene, due to low viscosity of SFCO<sub>2</sub>.

**Polymer Functionalization and Processing** – In 1994 and 1995, Watkins and McCarthy showed how SFCO<sub>2</sub> could be used to carry small molecules for functionalization of a polymer.


## Applications - CO<sub>2</sub> as reactant

In some reactions, CO<sub>2</sub> may be consumed, i.e. mitigated. For example, Noyori et al described the use of SFCO<sub>2</sub> reacting with hydrogen to make formic acid. It is possible to obtain higher reaction rates and longer catalyst lifetimes using this method.

## Implementation

### Case 1: Pressure Control


In this example, the pump is used in constant pressure mode. In constant pressure mode, the pump automatically displaces its volume to achieve the pressure requested by the operator. For SFCO<sub>2</sub>, pressure is directly related to solubility.



**Figure 2: Configuration for pressure control of SFCO<sub>2</sub> density and solubility during reaction**

### Case 2: Stable Mass Delivery

In this example, the pump is maintained at a constant temperature and is operated in constant flow mode. Steady back pressure is provided by the back pressure device. With known pressure, temperature, and volumetric displacement rate (“flow rate”), the mass of CO<sub>2</sub> delivered to the reactor is predictive.



**Figure 3: Configuration for stable mass delivery rate of CO<sub>2</sub> to reactor**

## Supercritical Fluid Pumps

For the applications described above, there are generally two types of pumps used: reciprocating and syringe. Reciprocating pumps have pistons with short strokes, so they need to refill frequently. Since fluid flow stops during refill, pressure fluctuations and density changes will result. This can cause unwanted precipitation of components, or other problems. Syringe pumps, considered pulseless, are better suited for these applications as the pressure and flow rates can be more accurately controlled. For most supercritical fluid applications, pressure must be maintained, so a constant pressure mode is needed. When pumps compress CO<sub>2</sub>,

even in the liquid state, heat generated will accumulate in the pump head. If this heat is not removed, incoming CO<sub>2</sub> could be inadvertently heated above the critical point, thereby impacting fill efficiencies. For proper operation, CO<sub>2</sub> pumps must incorporate some means to remove this heat.

Not all materials are suitable for use with CO<sub>2</sub>, so wetted materials must be checked for compatibility.

## Why Use Teledyne Isco Pumps?

Teledyne Isco syringe pumps are well suited for use with CO<sub>2</sub> and provide the best in accuracy and reliability. Flow rate accuracy is +/- 0.5% or better, and flows are pulseless. Pulseless flow means fluid pressure and density are constant, without changes in solvating properties. Pumps can be operated in either constant flow or constant pressure. For continuous operation, dual pump systems deliver fluid in unattended operation. Pumps can be operated as stand-alone or via external control. Teledyne Isco syringe pumps have a "poor fill alarm" which can alert the user if changing to a full supply bottle is needed.

Cooling jackets are available to maintain proper fluid temperature in the pump. Special valve packages for dual pump systems are CO<sub>2</sub> compatible.

## Recommendations for Teledyne Isco Pumps

Typically, chemical engineers who work with supercritical fluids choose to work with the Model 500D or 500HL pump. Sometimes, Model 260D or 100DM pumps are used in order to achieve higher pressure and/or more accurate flows at very slow flow rates. Single pumps are most often used in batch applications, while dual pumps are used in continuous flow.

**Table 1: Commonly Recommended Supercritical Fluid Pumps**

|                      | Model | Wetted Materials for CO <sub>2</sub> | Cooling/Heating Jacket | External Control Option |
|----------------------|-------|--------------------------------------|------------------------|-------------------------|
| Single Pump          | 260D  | Standard                             | Available              | Labview Driver          |
|                      | 500D  | Standard                             | Available              | Labview Driver          |
| Continuous Flow Pump | A260  | Standard                             | Available              | Labview Driver          |

**Table 2: Other Pump Models Available**

|                      | 1000D       | 500D        | 260D        | 100DX        | 100DM        | 65D          |
|----------------------|-------------|-------------|-------------|--------------|--------------|--------------|
| Flow Range (ml/min)  | 0.100 - 408 | 0.001 - 204 | 0.001 - 107 | 0.00001 - 60 | 0.00001 - 30 | 0.00001 - 25 |
| Pressure Range (psi) | 0 - 2,000   | 0 - 3,750   | 0 - 7,500   | 0 - 10,000   | 0 - 10,000   | 0 - 20,000   |

References:

- 1) Eric J. Beckman "Using CO<sub>2</sub> to Produce Chemical products Sustainably", *Environmental Science & Technology* (1 September 2002) 347A.
- 2) J. Calvin Giddings "High Pressure Gas Chromatography of Nonvolatile Species", *Science* 162 (4 October 1968) 67.
- 3) A.I. Cooper, et. al. "Extraction of a hydrophilic compound from water into liquid CO<sub>2</sub> using dendritic surfactants", *Nature* 389 (25 September 1997) 368.
- 4) P. A. Charpentier; J. M. DeSimone; G. W. Roberts "Continuous Polymerizations in Supercritical Carbon Dioxide." "In Clean Solvents," Abraham, M. A., Moens, L., Ed. ACS Symposium Series: Washington, DC, 2002; Vol. Chapter 9, pp 113-135.
- 5) J.B. McClain, et. al. "Design of Nonionic Surfactants for Supercritical Carbon Dioxide", *Science* 274 (20 December 1996) 2049.
- 6) Stephen K. Ritter "Microchips' Heavy Burden", *Chemical & Engineering News* (23 December 2002) 25.
- 7) David Filmore "The Greening of Catalysis", *Today's Chemist at Work* (November 2002) 29.
- 8) J.M. DeSimone, et. al. "Synthesis of Fluoropolymers in Supercritical Carbon Dioxide", *Science* 257 (14 August 1992) 945.
- 9) Jessop, P.G., Ikariya, T., Noyori, R. *Nature*, 1994, 368, 231.
- 10) Jessop, P.G., Hsiao, Y., Ikariya, T., Noyori, R. *J.Am.Chem. Soc.*, 199, 118, 344.
- 11) Liotta, C.L., Eckert, C.A., Hallett, J.P., Pollet, P. "Recycling Homogeneous Catalysts for Sustainable Technology", ORCS presentation, 2006.
- 12) Sparks, D.L., Hernandez, R., Zappi, M., French, T., Toghiani, H., Toghiani, R.K., Alley, E., "Oleic Acid Oxidation in Supercritical Carbon Dioxide", AIChE Spring presentation, 2006.
- 13) Sui, R; Rizkalla, A. S.; Charpentier, P. A. (2005) Formation of titania nanofibers: a direct sol-gel route in supercritical CO<sub>2</sub>. *Langmuir* 21 (14) 6150-6153.

*Last modified May 7, 2012*

**Teledyne Isco, Inc.**

P.O. Box 82531, Lincoln, Nebraska, 68501 USA  
Toll-free: (800) 228-4373 • Phone: (402) 464-0231 • Fax: (402) 465-3091  
E-mail: [IscoInfo@teledyne.com](mailto:IscoInfo@teledyne.com)

Teledyne Isco is continually improving its products and reserves the right to change product specifications, replacement parts, schematics, and instructions without notice.

