

CombiFlash[®] EZ Prep Installation Requirements

Company:

Shipping Address:	Contact:	Phone:
		Email:

Overview

This document provides general installation requirements and site preparation for the CombiFlash EZ Prep system.

Specifications

Table 1: Physical Specifications

CombiFlash EZ Prep		
H x W x D	68 x 43 x 43.2 cm	26.8 x 16.9 x 17"
Weight	43.6 kg	91.6 lb

Table 2: Laboratory Requirements

Environment	
Temperature	Recommended range: 20 - 30 °C Maximum range: 15 - 40 °C
Humidity	< 90%
Always maintain adequate ventilation to control vapors and ELSD exhaust.	
Electrical	
Nominal Voltage	± 10% (Example: 117VAC ± 11.7VAC)
One outlet is required.	
Solvents	
Four solvent capable.	
Gas Requirement- Nitrogen Recommended	
Purity	≥ 98%
Pressure	65± 5 psi (4.2-8.3 bar)
Gas Consumption	<2.51 l/min, max
System Connection	0.125" od, semi-rigid tubing
Space Requirement	
Ensure that adequate space is provided for the flash system solvent supply bottles and waste containers.	

Demo/Installation Supplies Required

- Nitrogen tank, house nitrogen line, or nitrogen generator

Note

The CombiFlash system needs to be less than 20 ft away from the nitrogen source.

- Two adjustable wrenches for nitrogen line hook up
- 65 psi gas pressure regulator connected to nitrogen
- Teflon tape
- 4 inlet solvents in 4 L bottles (Hexane, Ethyl Acetate, Water (HPLC Grade), Acetonitrile)
- Test tubes for rack (18x150 mm standard)
- Waste container- 4 L bottle
- Sample loop waste container- 100 ml flask
- Guard column and Prep column of choice
- 1 mL and 10 ml injection syringes for sample injection and sample loop wash
- 20 mL of IPA to fill P-trap on ELSD
- Secondary containment for solvent and waste containers, if necessary.

Created January 26, 2021

Use and Disclosure of Data: Information contained herein is classified as EAR99 under the U.S. Export Administration Regulations. Export, reexport or diversion contrary to U.S. law is prohibited.

Teledyne ISCO

P.O. Box 82531, Lincoln, Nebraska, 68501 USA
Toll-free: (800) 775-2965 • Phone: (402) 464-0231 • Fax: (402) 465-3001
E-mail: IscoService@teledyne.com

Teledyne ISCO is continually improving its products and reserves the right to change product specifications, replacement parts, schematics, and instructions without notice.

