

Improved Compound Recovery with ELSD Systems

Abstract

The ability to collect compounds is directly related to the ability to detect them. Many compounds possess weak or non-existent chromophores which causes UV detection to exhibit poor yields. The use of Evaporative Light Scattering Detection (ELSD) allows detection of these compounds; however, there are differences between the Teledyne Isco external ELSD systems and other systems. Compounds that are somewhat volatile can give a reduced response on the ELSD which reduces yield.

Background

Figure 1 below demonstrates how a weak signal reduces the amount of material collected because the material is not “visible” to the detector¹.





Figure 1: Shaded areas denote material recovered with an ELSD detector, but not detected with a UV detector

Although there are two peaks visible to the UV detector, the minor compound at ~5 column volumes (CV) is not detected by UV absorbance. The shaded areas on the ELSD peaks are also lost if only monitoring the UV signal. The sample is a mixture of tocopherols purified on a 30 g RediSep Rf Gold[®] diol column (PN 69-2203-516) with a hexane/ethyl acetate gradient. This example serves to demonstrate how compounds can be lost when a compound is only weakly detected.

Many flash ELSD systems exhibit poor response to semi-volatile compounds compared to the CombiFlash[®] Rf 200i (PN 68-5230-010) or the Model 340CF external ELSD (PN 68-5237-062); these systems exhibit reduced yield compared to the Teledyne Isco products.

Experimental and Results


A sample of 2,3-isopropylidene-D-ribofuranose was synthesized as previously described¹. The reaction mixture (1.0 mL) was purified on 40g RediSep[®] Rf silica columns (PN 69-2203-340) using the default gradient on each flash system; a hexanes/ethyl acetate gradient was run.



Fraction 1: 691.4 mg
 Fraction 2: 45.6 mg
 Fraction 3: 84.6 mg
Total recovery: 821.6mg

Figure 2: Model 340CF External ELSD

1. Bellinghausen, P; Fowler, N; Lewis, R; Silver, J. *Advanced Detection Techniques for Flash Chromatography*, Presented at 52nd Spring 2012 meeting of the American Chemical Society, San Diego, California, March 2012




Fraction 1: 692.1 mg

Fraction 2: 49.8 mg

Fraction 3: 37.6 mg

Total recovery: 821.6mg

Figure 3: CombiFlash Rf-200i with internal ELSD



Fraction 1: 444.6 mg

Total recovery: 444.6mg

Figure 4: Competitive flash system with ELSD

Conclusion

The CombiFlash systems demonstrated much greater recoveries than the competing systems- a yield improvement of nearly two-fold. The competing system detected only 64% of the major peak compared to the CombiFlash system, and failed to detect the minor components. For this example, the CombiFlash system provided 34% more compound for further synthesis.

Teledyne Isco

P.O. Box 82531, Lincoln, Nebraska, 68501 USA

Toll-free: (800) 228-4373 • Phone: (402) 464-0231 • Fax: (402) 465-3091

E-mail: IscoInfo@teledyne.com

Teledyne Isco is continually improving its products and reserves the right to change product specifications, replacement parts, schematics, and instructions without notice.

