

RediSep® Rf Cartridges

NEW

RediSep® Rf

easy-load cartridges

- Convenient volume to introduce your compound
- Better for solid samples
- Doesn't dilute the sample introduction
- Reduced band broadening for smaller samples
- Compatible with all Teledyne Isco sample load cartridge caps

Cartridge Size	Maximum Suggested Sample Load
2.5g	0.25 - 0.5g
5g	0.5 - 1g
12g	1 - 2.5g
25g	2.5 - 5g
32g	5 - 6.5g
65g	6.5 - 13g
12.5g	13 - 25g
260g	25 - 72g



Sample Load Prepacked Silica Gel Rf Cartridges

2.5-gram, pkg. of 20	69-3873-308	NEW
(For use on all Combiflash systems except OptiX 10.)		
5-gram, pkg. of 20	69-3873-219	
(For use on all Combiflash systems except OptiX 10, see 69-3873-145.)		
12-gram, pkg. of 15	69-3873-309	NEW
25-gram, pkg. of 15	69-3873-221	
32-gram, pkg. of 12	69-3873-310	NEW
65-gram, pkg. of 12	69-3873-226	

Sample Load Prepacked Silica Gel Cartridges for Companion XL

125-gram, pkg. of 4	69-3873-311	NEW
260-gram, pkg. of 4	69-3873-202	

Sample Load Prepacked Silica Gel Cartridges for OptiX 10

5-gram, pkg. of 20	69-3873-145
(For use with OptiX 10 systems only. Use 69-3873-219 for other systems.)	