


# RediSep Rf Gold™ — High Performance Gram - Multigram Purification

## Save Time & Solvent

The CombiFlash® Torrent™ system and new RediSep Rf Gold™ High Performance silica columns team together to bring Flash chromatography to the next level. Take advantage of the sharper peaks provided by spherical media to shorten purification time. Convert your methods to Gold Speed with shorter gradients.

- ◆ Save hundreds of dollars with half the solvent usage
- ◆ Reduce purification time up to 40%
- ◆ Reduce flammable solvent requirements
- ◆ Separate silica sensitive compounds faster



For additional information, contact us today! **(800) 228-4373**



## RediSep Rf vs. RediSep Rf Gold

Column Size	Sample Size	Solvent Usage	
		Gold Resolution	Gold Speed
80 g	0.8 - 8 g	2.4 L	1.4 L
120 g	1.2 - 12 g	3.7 L	2.2 L
220 g	2.2 - 22 g	6.5 L	3.8 L
330 g	3.3 - 33 g	9 L	5 L
750 g	7.5 - 75 g	23 L	13 L
1.5 kg	15 - 150 g	46 L	26 L
3 kg	30 - 300 g	92 L	53 L

## Ordering Information

### RediSep Rf Gold™ Silica Columns

120 gram, package of 6 .....	.69-2203-349
220 gram, package of 4 .....	.69-2203-359
330 gram, package of 3 .....	.69-2203-369
750 gram, package of 3 .....	.69-2203-427
1.5 kilogram, package of 2 .....	.69-2203-428
3 kilogram, package of 1 .....	.69-2203-529

These RediSep Rf Gold™ columns are packed with 20-40 µm, spherical gel.

U.S. patent #7,008,541 and European patent #1,316,798 granted.

Teledyne Isco, Inc. reserves the right to change specifications without notice.

To learn more about RediSep Rf Gold™ columns, visit our  
Liquid Chromatography Application Note Library at  
[www.isco.com/lcappnotes](http://www.isco.com/lcappnotes).



**TELEDYNE ISCO**

A Teledyne Technologies Company

4700 Superior Street • Lincoln, NE 68504 USA

Toll free: (800) 228-4373 (USA & Canada) • Phone: (402) 464-0231 • Fax: (402) 465-3022

e-mail: [iscoinfo@teledyne.com](mailto:iscoinfo@teledyne.com)