

ReaXus[®] M1 Class

Single Piston Pumps with Excellent Capability at a Low Cost



The compact M1 Class consists of single-headed, positive displacement piston pumps that deliver high-performance and precision at an affordable price. Superior to Gear, Diaphragm, and Peristaltic pumps, the M1 Class provides

Superior to Gear, Diaphragm, and Peristaltic pumps, the M1 Class provides greater performance at an equivalent cost to lesser-functioning products.

greater performance at an equivalent cost to lesser-functioning products.



Applications include precise chemical metering, dispensing, and specific HPLC

separations. Standard fluid path materials are Stainless Steel and PEEK. With 10 mL/min, 40 mL/min, and 100 mL/min versions, reaching pressures up to 2,500 psi, these pumps capture the majority of typical laboratory fluid transfer applications.

Features include interactive keypad control, Outlet Filter, RS-232 Serial and Micro USB 2.0 ports for complete PC control and status, as well as hardware for laboratory stand mounting. Having the smallest possible footprint for a high pressure pump, the M1 Class will easily integrate into any laboratory or industrial environment.

REAXUS[®]

FLOW RATE

10 mL/min
40 mL/min
100 mL/min

PRESSURE

Up To 2,500 psi
(10 mL/min)

FLUID PATH

Stainless Steel or PEEK

CONTROL

RS-232, Micro USB, Run/Stop

M1 Class Specifications

	10 mL/min	40 mL/min	100 mL/min
Flow Rate Range	0.01 – 10.000 mL/min	0.1 – 40.0 mL/min	0.1 – 100.0 mL/min
Flow Rate Increment	0.01 mL/min	0.1 mL/min	0.1 mL/min
Flow Rate Accuracy*	±5% (0.5 – 10.0 mL/min; 80:20 Water/IPA; 1,000 psi)	±5% (1.0 – 36.0 mL/min; 80:20 Water/IPA; 250 psi)	±5% (10.0 – 100.0 mL/min; 80:20 Water/IPA; 100 psi)
Flow Reproducibility	1.5% RSD (0.3 – 10.0 mL/min; 80:20 Water/IPA; 1,000 psi)	1% RSD (1.0 – 36.0 mL/min; 80:20 Water/IPA; 250 psi)	1% RSD (2.0 – 100.0 mL/min; 80:20 Water/IPA; 100 psi)
Maximum Delivery Pressure	2,500 psi	500 psi	250 psi
Pressure Accuracy	N/A		
Piston Displacement	25.1 µL	100 µL	226 µL
Piston Wash	Self-flushing pump heads – continuous wash without auxiliary pump		N/A
Wetted Materials	PEEK pumps may include: PEEK, Synthetic Ruby, Sapphire, UHMWPE, FFKM Stainless Steel pumps may include: 316 Stainless Steel, PEEK, Synthetic Ruby, Sapphire, Alumina, UHMWPE, PTFE, ETFE		
Dimensions (H x W x D)	5.6" x 3.0" x 10.6" (14.2 x 7.6 x 26.9 cm)		
Weight	6 lbs (2.7 kg)		
Power	100 – 240 Vac (±10%), 50 – 60 Hz, 40 W		
Front Panel	4-digit LED to monitor flow rate and indicate alarms. Run/Stop and menu keys, two status indicators, Increase/Decrease displayed parameters.		
Communication	RS-232, micro-USB, Run/Stop TTL		
Viscosity**	0.5 – 7.0 cP	0.5 – 5.2 cP	0.5 – 6.0 cP
Liquid Temperature	0 – 82 °C		
Environment	Indoor use only (Pollution degree 2)		
Altitude	Up to 2000 m (6500 ft)		
Ambient Temperature	10 – 30 °C		
Humidity	20 – 90% relative humidity		

*Flow Rate is dependent on solvent selection and operating pressure

**Range of fluid viscosities where published flow accuracy is maintained

M1 CLASS—Order Numbers

CONSTANT FLOW

STAINLESS STEEL

Stainless Steel Fluid Path; with Pulse Dampener; Standard Seal (Stainless Steel Energizer); Constant-Flow Control Mode

0.01–10.00 mL/min; 2,000 psi # M1010SFP1XFTL

0.01–40.0 mL/min; 500 psi # M1040SFP1XFTL

0.01–100.0 mL/min; 100 psi # M1100SNN1XFTL

Teledyne ISCO

P.O. Box 82531, Lincoln, Nebraska, 68501 USA
Toll-free: (800) 228-4373 • Phone: (402) 464-0231 • Fax: (402) 465-3091

teledyneisco.com



Teledyne ISCO is continually improving its products and reserves the right to change product specifications, replacement parts, schematics, and instructions without notice.



L-7149 Rev 7.0
08/23